



Recombinant *Saccharomyces cerevisiae* General negative regulator of transcription subunit 2 (CDC36)

Product Code	CSB-MP361710SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06100
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MEKFGKALV PLLKLEDKEL SSTYDHSMTL GADLSSMLYS LGIPRDSQDH RVLDTFQSPW AETSRSEVEP RFFTPESFTN IPGVLQSTVT PPCFNSIQND QQRVALFQDE TLFFLFYKHP GTVIQELTYL ELRKRNWRYH KTLKAWLTKD PMMEPIVSAD GLSERGSYVF FDPQRWEKCQ RDFLLFYNAI M
Source	Mammalian cell
Target Names	CDC36
Protein Names	Recommended name: General negative regulator of transcription subunit 2 Alternative name(s): cell division cycle protein 36
Expression Region	1-191
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.