



Recombinant Solanum tuberosum Photosystem II 10 kDa polypeptide, chloroplastic (PSBR)

Product Code	CSB-BP361738FIG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06183
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	S GVKKCLKTDKP YGINGSMALR DGVDASGRKP KGKGVYQYVD KYGANVDGYS PIYNTDEWSP SGDVYVGGTT GLAIWAVTLV GILAGGALLV YNTSALAQ
Source	Baculovirus
Target Names	PSBR
Protein Names	Recommended name: Photosystem II 10 kDa polypeptide, chloroplastic Alternative name(s): Light-inducible tissue-specific ST-LS1 protein
Expression Region	40-138
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.