



Recombinant *Saccharomyces cerevisiae* cAMP-dependent protein kinase type 1 (TPK1)

Product Code	CSB-EP361759SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06244
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	STEEQNGGG QKSLDDRQGE ESQKGETSER ETTATESGNE SKSVEKEGGE TQEKPQPHV TYYNEEQYKQ FIAQARVTSG KYSLQDFQIL RTLGTGSFGR VHLIRSRHNG RYYAMKVLKK EIVVRLKQVE HTNDERLMLS IVTHPFIIRM WGTFQDAQQI FMIMDYIEGG ELFSLLRKSQ RFPNPVAKFY AAEVCLALEY LHSDIIRYRD LKPENILLDK NGHIKITDFG FAKYVPDVTY TLCGTPDYIA PEVVSTKPYN KSIDWWSFGI LIYEMLAGYT PFYDSNTMKT YEKILNAELR FPPFFNEDVK DLLSRLITRD LSQRLGNLQN GTEDVKNHPW FKEVVWEKLL SRNIETPYEP PIQQGQGDTS QFDKYPEEDI NYGVQGEDPY ADLFRDF
Source	E.coli
Target Names	TPK1
Protein Names	Recommended name: cAMP-dependent protein kinase type 1 Short name= PKA 1 EC= 2.7.11.11 Alternative name(s): CDC25-suppressing protein kinase PK-25
Expression Region	2-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.