



# Recombinant Human C-X-C motif chemokine 14 (CXCL14)

<b>Product Code</b>	CSB-BP006244HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O95715
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SKCKCS RKGPKIRYSV VKKLEMKPKY PHCEEKVMVII TTKSVSRVYRG QEHCLHPKLQ STKRFIKWYN AWNEKRRVYE E
<b>Source</b>	Baculovirus
<b>Target Names</b>	CXCL14
<b>Protein Names</b>	Recommended name: C-X-C motif chemokine 14 Alternative name(s): Chemokine BRAK MIP-2G Small-inducible cytokine B14
<b>Expression Region</b>	35-111
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene belongs to the cytokine gene family which encode secreted proteins involved in immunoregulatory and inflammatory processes. This protein is structurally related to the CXC (Cys-X-Cys) subfamily of cytokines. Members of this subfamily are characterized by two cysteines separated by a single amino acid. This cytokine displays chemotactic activity for monocytes but not for lymphocytes, dendritic cells, neutrophils or macrophages. It has been implicated that this cytokine is involved in the homeostasis of monocyte-derived macrophages rather than in inflammation.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.