



Recombinant Human Polypyrimidine tract-binding protein 3 (PTBP3)

Product Code	CSB-EP020060HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95758
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDGVVTDLIT VGLKRGSEDEL LSSGIINGPF TMNSSTPSTA NGNDSKKFKR DRPPCSPSRV LHLRKIPCDV TEAEIISLGL PFGKVTNLLM LKGKSQAFLE MASEEAAVVM VNYYPITPH LRSQPVYIQY SNHRELKTDN LPNQARAQAA LQAVSAVQSG SLALSGGPSN EGTVLPQGSP VLRIIIENLF YPVTLEVLHQ IFSKFGTVLK IITFTKNNQF QALLQYADPV NAHYAKMALD GQNIYNACCT LRIDFSKLTSLNVKYNNDKSRDFTRLDLPT GDGQPSLEPP MAAAFGAPGI ISSPYAGAAG FAPAIGFPQA TGLSVPVAVPG ALGPLTITSS AVTGRMAIPG ASGIPGNSVL LVTNLNPDLI TPHGLFILFG VYGDVHRVKI MFNKKENALV QMADANQAQL AMNHLGQRL YGKVLRLATLS KHQAVQLPRE GQEDQGLTKD FSNSPLHRFK KPGSKNFQNI FPPSATLHLS NIPPSVTVD LKNLFI EAGC SVKAFKFFQK DRKMALIQLG SVEEAIQALI ELHNHDLGEN HHLRVSFSSKSTI
Source	E.coli
Target Names	PTBP3
Protein Names	Recommended name: Polypyrimidine tract-binding protein 3 Alternative name(s): Regulator of differentiation 1 Short name= Rod1
Expression Region	1-552
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.