



# Recombinant Human Polypyrimidine tract-binding protein 3 (PTBP3)

<b>Product Code</b>	CSB-MP020060HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O95758
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDGVVTDLIT VGLKRGSEDEL LSSGIINGPF TMNSSTPSTA NGNDSKKFKR DRPPCSPSRV LHLRKIPCDV TEAEIISLGL PFGKVTNLLM LKGKSQAFLE MASEEAAVTM VNYYPITPH LRSQPVYIQY SNHRELKTDN LPNQARAQAA LQAVSAVQSG SLALSGGPSN EGTVLPGQSP VLRIIIENLF YPVTLEVLHQ IFSKFGTVLK IITFTKNNQF QALLQYADPV NAHYAKMALD GQNIYNACCT LRIDFSKLTLS LNVKYNNDKS RDFTRLDLPT GDGQPSLEPP MAAAFGAPGI ISSPYAGAAG FAPAIGFPQA TGLSVPVAVPG ALGPLTITSS AVTGRMAIPG ASGIPGNSVL LVTNLNPDLI TPHGLFILFG VYGDVHRVKI MFNKKENALV QMADANQAQL AMNHLGQRL YGKVLRLATLS KHQAVQLPRE GQEDQGLTKD FSNSPLHRFK KPGSKNFQNI FPPSATLHLS NIPPSVTVDD LKNLFI EAGC SVKAFKFFQK DRKMALIQLG SVEEAIQALI ELHNHDLGEN HHLRVSFSKS TI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PTBP3
<b>Protein Names</b>	Recommended name: Polypyrimidine tract-binding protein 3 Alternative name(s): Regulator of differentiation 1 Short name= Rod1
<b>Expression Region</b>	1-552
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.