



Recombinant *Drosophila melanogaster* Tropomyosin-1, isoforms 9A/A/B (Tm1)

Product Code	CSB-MP361937DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06754
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MDAIKKKMQA MKVDKDGAL RALVCEQEAR DANTRA EKAE EEARQLQKKI QTVENELDQT QEALTLVTGK LEEKNKALQN KKKTTKMTTS IPQGTL LDVL KKKMRQTKEE MEKYKDECEE FHKRLQLEVV RREEAESEVA ALNRRIQ LLE EDLERSEERL GSATAKLSEA SQA ADESERA RKILEN RALA DEERMDALEN QLKEARFLAE EADKKYDEVA RKLAMVEADL ERAEERAEQG ENKIVELEEE LRVVGNNLKS LEVSEEKANQ REEYKNQIK TLNTRLKEAE ARAEFAERSV QKLQKEVDRL EDDLVL EKER YKDIGDDLDT AFVELILKE
Source	Mammalian cell
Target Names	Tm1
Protein Names	Recommended name: Tropomyosin-1, isoforms 9A/A/B Alternative name(s): Cytoskeletal tropomyosin Tropomyosin II
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.