



Recombinant *Saccharomyces cerevisiae* Protein SPT3 (SPT3)

Product Code	CSB-MP361965SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06844
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MMDKHKYRVE IQQMMFVSGE INDPPVETTS LIEDIVRGQV IEILLQSNKT AHLRGSRSIL PEDVIFLIRH DKAKVNRLRT YLSWKDLRKN AKDQDASAGV ASGTGNPGAG GEDDLKKAGG GEKDEKDGGN MMKVKKSQIK LPWELQFMFN EHPLENNDDN DDMEDEDEREA NIVTLKRLKM ADDRTRNMTK EEYVHWSDCR QASFTFRKNK RFKDWSGISQ LTEGKPHDDV IDILGFLTFE IVCSLTETAL KIKQREQVLQ TQKDKSQQSS QDNTNFEFAS STLHRKKRLF DGPENVINPL KPRHIEEAWR VLQTIDMRHR ALTNFKGGRL SSKPIIM
Source	Mammalian cell
Target Names	SPT3
Protein Names	Recommended name: Protein SPT3 Alternative name(s): Positive regulator of Ty transcription
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.