



Recombinant Gag polyprotein (ev-1)

Product Code	CSB-YP361990ARE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06936
Product Type	Recombinant Protein
Immunogen Species	Avian endogenous virus EV-1
Purity	>85% (SDS-PAGE)
Sequence	MEAVIKVISS ACKTYCGKTS PSKKEIGAML SLLQKEGLLM SPSDLYSPGS WDPITAALSQ RAMVLGKSGE LKTWGLVLGA LKAAREEQVT SEQAKFWLGL GGGRVSPPGP ECIEKPATER RIDKGEEVGE TTAQRDAKMA PEKMATPKTV GTSCYQCGTA TGCNCATASA PPPPYVG
Source	Yeast
Target Names	ev-1
Protein Names	Recommended name: Gag polyprotein Alternative name(s): Core polyprotein Cleaved into the following 4 chains: 1. Matrix protein p19 2. Core protein p23 3. Core protein p2 4. Core protein p10
Expression Region	1-177
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.