



Recombinant Enterobacteria phage P1 Recombinase cre (cre)

Product Code	CSB-MP362000EDD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06956
Product Type	Recombinant Protein
Immunogen Species	Escherichia phage P1 (Bacteriophage P1)
Purity	>85% (SDS-PAGE)
Sequence	MSNLLTVHQN LPALPVDATS DEVRKNLMDM FRDRQAFSEH TWKMLLSVCR SWAAWCKLNN RKFPAEPED VRDYLLYLQA RGLAVKTIQQ HLGQLNMLHR RSGLPRPSDS NAVSLVMRRI RKENVDAGER AKQALAFERT DFDQVRSIME NSDRCQDIRN LAFLGIAYNT LLRIAEIARI RVKDISRTDG GRMLIHIGRT KTLVSTAGVE KALSLGVTKL VERWISVSGV ADDPNNYLFC RVRKNGVAAP SATSQLSTRA LEGIFEATHR LIYGAKDDSG QRYLAWSGHS ARVGAARDMA RAGVSIPEIM QAGGWTVNI VMNYIRNLDS ETGAMVRLLE DGD
Source	Mammalian cell
Target Names	cre
Protein Names	Recommended name: Recombinase cre
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.