



Recombinant Escherichia coli 6-phosphofructokinase isozyme 2 (pfkB)

Product Code	CSB-MP362011ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06999
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MVRIYTLTLA PSLDSATITP QIYPEGKLR C TAPVFEPGGG GINVARIAIH LGG SATAIFP AGGATGEHLV SLLADENVPV ATVEAKDWTR QNLHVHVEAS GEQYRFVMPG AALNEDEF RQ LEEQVLEIES GAILVISGSL PPGVKLEKLT QLISAAQKQG IRCIVDSSGE ALSAAL AIGN IELVKPNQKE LSALVNRELT QPDDVRKAAQ EIVNSGKAKR VV VSLGPQGA LGVDS ENCIQ VVPPP VKS QS TVGAGDSMVG AMTLKLAENA SLEEMVRFGV AAGSAATLNQ GTRLC SHDDT QKIYAYLSR
Source	Mammalian cell
Target Names	pfkB
Protein Names	Recommended name: 6-phosphofructokinase isozyme 2 Short name= Phosphofructokinase-2 EC= 2.7.1.11
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.