



Recombinant Enterobacteria phage T4 DNA topoisomerase medium subunit (52)

Product Code	CSB-MP362030EDZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07065
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	>85% (SDS-PAGE)
Sequence	MQLNNRDLKS IIDNEALAYA MYTVENRAIP NMIDGFKPVQ RFVIARALDL ARGNKDKFHK LASIAGGVAD LGYHHGENSA QDAGALMANT WNNNFPLLDG QGNFGSRTVQ KAAASRYIFA RVSKNFYINVY KDTEYAPVHQ DKEHIPPAFY LPIIPTVLLN GVSGIATGYA TYILPHSVSS VKKAVLQALQ GKKVTKPKVE FPEFRGEVVE IDGQYEIRGT YKFTSRTQMH ITEIPYKYDR ETYVSKILD LENKGFITWD DACGEHGF GF KVKFRKEYSL SDNEEERHAK IMKDFGLIER RSQNITVINE KGKLQVYD NV VDLIKDFVEV RKTYYQKRID NKIKETESAF RLAFAKAHFI KKVISGEIVV QGKTRKELTE ELSKIDMYSS YVDKLVGMNI FHMTSDEAKK LAEEAKAKKE ENEYWKTTDV VTEYTKDLEE IK
Source	Mammalian cell
Target Names	52
Protein Names	Recommended name: DNA topoisomerase medium subunit EC= 5.99.1.3 Alternative name(s): Protein Gp52
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.