




Recombinant Human rhinovirus A serotype 89 Genome polyprotein, partial

Product Code	CSB-EP362073HQD
Abbreviation	Recombinant Human rhinovirus A serotype 89 Genome polyprotein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07210
Product Type	Recombinant Protein
Immunogen Species	Human rhinovirus A serotype 89 (strain 41467-Gallo) (HRV-89)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	NPVENYIDSVLNEVLVVPNIQPSTSVSSHAAPALDAAETGHTSSVQPEDMIETR YVITDQTRDETSIESFLGRSGCIAMIEFNTSSDKTEHDKIGKGFKTWKVSLQEM AQIRRKYEFLTYTRFDSEITIVTAAAQGNDSGHIVLQFMYVPPGAPVPEKRDD YTWQSGTNASVFWQEGQPYPRTIPFMSIASAYYMFYDGYDGDSSAASKYGSV VTNDMGTCVRIVTSNQKHDSNIVCRIYHKAKHIKAWCPRPPRAVAYQHTHSTN YIPSNGEATTQIKTRPDVFTVTNV
Source	E.coli
Expression Region	575-866aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	48.6 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.