



Recombinant Escherichia coli Chemotaxis protein CheA (cheA)

Product Code	CSB-MP362114ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07363
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSMDISDFYQ TFFDEADELL ADMEQHLLVL QPEAPDAEQL NAIFRAAHSI KGGAGTFGFS VLQETTHLME NLLDEARRGE MQLNTDIINL FLETKDIMQE QLDAYKQSQE PDAASFDYIC QALRQLALEA KGETPSAVTR LSVVAKSEPO DEQSRSQSPR RIILSRLKAG EVDLLEEELG HLTTLTDVVK GADSLSAILP GDIAEDDITA VLCFVIEADQ ITFETVEVSP KISTPPVLKL AAEQAPTGRV EREKTTTRSNE STSIRVAVEK VDQLINLVGE LVITQSMLAQ RSSELDPVNH GDLITSMGQL QRNARDLQES VMSIRMM PME YVFSRYPRLV RDLAGKLGKQ VELTLVGSST ELDKSLIERI IDPLTHLVRN SLDHGIELPE KRLAAGKNSV GNLILSAEHQ GGNICIEVTD DGAGLNRERI LAKAASQGLT VSENMSDDEV AMLIFAPGFS TAEQVTDVSG RGVGMDVVKR NIQKMGGHVE IQSKQGTGTT IRILLPLTLA ILDGMSVRVA DEVFILPLNA VMESLQPREA DLHPLAGGER VLEVRGEYLP IVELWKVFNV AGAKTEATQG IVVILQSGGR RYALLVDQLI GQHQVVVKNL ESNYRKVPGI SAATILGDGS VALIVDVSAL QAINREQRMA NTAA
Source	Mammalian cell
Target Names	cheA
Protein Names	Recommended name: Chemotaxis protein CheA EC= 2.7.13.3
Expression Region	1-654
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.