



Recombinant *Saccharomyces cerevisiae* Ras-related protein SEC4 (SEC4)

Product Code	CSB-YP362164SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07560
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSGLRTVSAS SGNGKSYDSI MKILLIGDSG VGKSCLLVRF VEDKFNPSFI TTIGIDFKIK TVDINGKKVK LQLWDTAGQE RFRTITTAYY RGAMGIILVY DVTDERFTFN IKQWFKTVNE HANDEAQLLL VGNKSDMETR VVTADQGEAL AKELGIPFIE SSAKNDDNVN EIFFTLAKLI QEKIDSNKLV GVGNGKEGNI SINSGSGNSS KSNCC
Source	Yeast
Target Names	SEC4
Protein Names	Recommended name: Ras-related protein SEC4 Alternative name(s): Suppressor of RHO3 protein 6
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.