



Recombinant Plasmid partition protein A (parA)

Product Code	CSB-EP362181ENL-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07620
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	MSDSSQLHKV AQRANRMLNV LTEQVQLQKD ELHANEFYQV YAKAALAKLP LLTRANVDYA VSEMEEKGYV FDKRPAGSSM KYAMSIQNI DIYEHARGVPK YRDRYSEAYV IFISNLKGGV SKTVSTVSLA HAMRAHPHLL MEDLRILVID LDPQSSATMF LSHKHSIGIV NATSAQAMLQ NVSREELLEE FIVPSVVPGV DVMPASIDDA FIASDWRELC NEHLPQQNIH AVLKENVIDK LKSDYDFILV DSGPHLDAFL KNALASANIL FTPLPPATVD FHSSLKYVAR LPELVKLISD EGCECQLATN IGFMSKLSNK ADHKYCHSLA KEVFGGDMLD VFLPRLDGFE RCGESFDTVI SANPATYVGS ADALKNARIA AEDFAKAVFD RIEFIRSN
Source	E.coli
Target Names	parA
Protein Names	Recommended name: Plasmid partition protein A
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.