



Recombinant Escherichia coli Homoserine O-succinyltransferase (metA)

Product Code	CSB-BP362182ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07623
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	PIRVPDELP AVNFLREENV FVMTTSRASG QEIRPLKVL I LNLMPKKIET ENQFLRLLSN SPLQVDIQLL RIDSRESRNT PAEHLNNFYC NFEDIQDQNF DGLIVTGAPL GLVEFNDVAY WPQIKQVLEW SKDHVTSTLF VCWAVQAALN ILYGIPKQTR TEKLSGVYEH HILHPHALLT RGFDDSFAP HSRYADFPAA LIRDYTDLEI LAETEEGDAY LFASKDKRIA FVTGHPEYDA QTLAQEFFRD VEAGLDPDVP YNYFPHNDPQ NTPRASWRSH GNLLFTNWLN YVYQITPYD LRHMNPTLD
Source	Baculovirus
Target Names	metAS
Protein Names	Recommended name: Homoserine O-succinyltransferase EC= 2.3.1.46 Alternative name(s): Homoserine O-transsuccinylase Short name= HTS
Expression Region	2-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.