



Recombinant Enterobacteria phage K3 Tail fiber protein p36 (36)

Product Code	CSB-EP362183ECV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07626
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage K3 (Bacteriophage K3)
Purity	>85% (SDS-PAGE)
Sequence	MADLKVGGSTT GGSVIWHQGN FPLNPAGDDV LYKSFKIYSE YNKPQATDND FVSKANGGTK SQGFQKEVEF REGVKISATF SGGSDLNGLY SGNGDGASRE KANMDLRSWY GIGIWNTCTG DGSGTHNSGL MPGVVHVETC WRLLFQTDNV ISNASGPTGP AHLTRKDYVD GAINVTANA NSRVLRSGDT MTGNLTAPNF FSQNPASQPS HVPRFDQIVI KDSVQDFGY
Source	E.coli
Target Names	36
Protein Names	Recommended name: Tail fiber protein p36 Alternative name(s): Protein Gp36
Expression Region	1-230
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.