



Recombinant *Drosophila melanogaster* Serendipity locus protein delta (Sry-delta)

Product Code	CSB-MP362191DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07664
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MDTCFFCGAV DLSDTGSSSS MRYETLSAKV PSSQKTVSLV LTHLANCIQT QLDLKPGARL CPRCFQELSD YDTIMVNLMT TQKRLTTQLK GALKSEFEVP ESGEDILVEE VEIPQSDVET DADAEADALF VELVKDQEEES DTEIKREFVD EEEEEDDDDD DEFICEDVDV GDSEALYGKS SDGEDRPTKK RVKQECTTCG KVYNSWYQLQ KHISEEHSKQ PNHICPICGV IRRDEEYLEL HMNLHEGKTE KQCRYCPKSF SRPVNTRLRHM RMHWDKKKYQ CEKCGLRFSQ DNLLYNHRLR HEAEENPIIC SICNVSFKSR KTFNHHTLIH KENRPRHYCS VCPKSFTERY TLKMHMKTHE GDVVYGVREE APADEQQVVE ELHVDVDESE AAVTVIMSDN DENSGFCLIC NTNFNKEL EHLQFDHDV VLK
Source	Mammalian cell
Target Names	Sry-delta
Protein Names	Recommended name: Serendipity locus protein delta
Expression Region	1-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.