



Recombinant Schizosaccharomyces pombe 60S acidic ribosomal protein P2-alpha (rpp201)

Product Code	CSB-MP362322SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08094
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKYLAAYLLL TVGGKDSPSA SDIESVLSTV GIEAESERIE TLINELNGKD IDELIAAGNE KLATVPTGGA ASAAPAAAAG GAAPAAEEEA KEEAKEEEEE DEDMGFGLFD
Source	Mammalian cell
Target Names	rpp201
Protein Names	Recommended name: 60S acidic ribosomal protein P2-alpha Short name= A2 Alternative name(s): L12EI L40C
Expression Region	1-110
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.