



Recombinant *Nicotiana tabacum* Endochitinase A (CHN48)

Product Code	CSB-BP362364NHE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08252
Product Type	Recombinant Protein
Immunogen Species	<i>Nicotiana tabacum</i> (Common tobacco)
Purity	>85% (SDS-PAGE)
Sequence	EQCGSQA GGARCPSGLC CSKFGWCGNT NDYCGPGNCQ SQCPGGPTPT PPTPPGGGDL GSISSSMFD QMLKHRNDNA CQGKGFYSYN AFINAARSFP GFGTSGDTTA RKREIAAFFA QTSHETTGGW ATAPDGPYAW GYCWLREQGS PGDYCTPSGQ WPCAPGRKYF GRGPIQISHN YNYGPCGRAI GVDLLNPNPD VATDPVISFK SALWFWMTPQ SPKPSCHDVI IGRWQPSAGD RAANRLPGFG VITNIINGGL ECGRGTDSDRV QDRIGFYRRY CSILGVSPGD NLDCGNQRSF GN
Source	Baculovirus
Target Names	CHN48
Protein Names	Recommended name: Endochitinase A Short name= CHN-A EC= 3.2.1.14
Expression Region	24-322
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.