



# Recombinant Human Zinc finger protein 264 (ZNF264)

<b>Product Code</b>	CSB-EP026622HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O43296
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAAAVLTDRA QVSVTFDDVA VTFTKEEWGQ LDLAQRTLYQ EVMLENCGLL VSLGCPVPKA ELICHLHGQ EPWTRKEDLS QDTCPGDKGK PKTTEPTTCE PALSEGISLQ GQVTQGNSVD SQLGQAEDQD GLSEMQEGHF RPGIDPQEK PGKMSPECDG LGTADGVCSR IGQEQVSPGD RVRSHNSCES GKDPMIQEEE NNFKCSECGK VFNKKHLLAG HEKIHSGVVP YECTECGKTF IKSTHLLQHH MIHTGERPYE CMECGKAFNR KSYLTQHQR I HSGEKPYKCN ECGKAFTHRS NFVLHNRHT GEKSFVCTEC GQVFRHRPGF LRHYVVHSGE NPYECLECGK VFKHRSYLMW HQTHTGEKP YECSECGKVF LESAALIHYY VIHTGEKPF CLECGKAFNH RSYLKRHQRI HTGEKPFVCS ECGKAFTHCS TFILHKRAHT GEKPFECCK GKAFSNRKDL IRHFSIHTGE KPYECVECGK AFTRMSG LTR HKRIHSGEKP YECVECGKSF CWSTNLIRHA IIHTGEKPYK CSECGKAFSR SSSLTQHQR HTGKNPISVT DVGRPFTSGQ TSVTLRELLL GKDFLNV TTE ANILPEETSS SASDQPYQRE TPQVSSL
<b>Source</b>	E.coli
<b>Target Names</b>	ZNF264
<b>Protein Names</b>	Recommended name: Zinc finger protein 264
<b>Expression Region</b>	1-627
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.