



Recombinant Human Zinc finger protein 264 (ZNF264)

Product Code	CSB-YP026622HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43296
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MAAAVLTDRA QVSVTFDDVA VTFTKEEWGQ LDLAQRTLYQ EVMLENCGLL VSLGCPVPKA ELICHLHGQ EPWTRKEDLS QDTCPGDKGK PKTTEPTTCE PALSEGISLQ GQVTQGNSVD SQLGQAEDQD GLSEMQEGHF RPGIDPQEK PGKMSPECDG LGTADGVCSR IGQEQVSPGD RVRSHNSCES GKDPMIQEEE NNFKCSECGK VFNKKHLLAG HEKIHSGVVP YECTECGKTF IKSTHLLQHH MIHTGERPYE CMECGKAFNR KSYLTQHQR IHSGEKPYKCN ECGKAFTHS NFVLHNRHT GEKSFVCTEC GQVFRHRPGF LRHYVVHSGE NPYECLECGK VFKHRSYLMW HQTHTGEKP YECSECGKVF LESAALIHYY VIHTGEKPF CLECGKAFNH RSYLKRHQRI HTGEKPFVCS ECGKAFTHCS TFILHKRAHT GEKPFECCK GKAFSNRKDL IRHFSIHTGE KPYECVECGK AFTRMSG LTR HKRIHSGEKP YECVECGKSF CWSTNLIRHA IIHTGEKPYK CSECGKAFSR SSSLTQHQR HTGKNPISVT DVGRPFTSGQ TSVTLRELLL GKDFLNV TTE ANILPEETSS SASDQPYQRE TPQVSSL</p>
Source	Yeast
Target Names	ZNF264
Protein Names	Recommended name: Zinc finger protein 264
Expression Region	1-627
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.