



Recombinant *Saccharomyces cerevisiae* Mitochondrial nuclease (NUC1)

Product Code	CSB-EP362435SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08466
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MCSRILLSGL VGLGAGTGLT YLLLNKHSPT QIIETPYPT QKPNSNIQSH SFNVDPGFF KYGFGPIHD LQNREEFISC YNRQTQNPYW VLEHITPESL AARNADRKNS FFKEDEVIPE KFRGKLRDYF RSGYDRGHQA PAADAKFSQQ AMDDTFYLSN MCPQVGEGFN RDYWAHLEYF CRGLTKKYKS VRIVTGPLYL PKKDPIDNKF RVNYEVIGNP PSIAVPTHFF KLIVAEAPTA NPAREDIAVA AFVLPNEPIS NETKLDFEV PIDALERSTG LELLQKVPPS KKKALCKEVN CQIVVRDFSN AAIKQSKDVK LLPPPCKRN
Source	E.coli
Target Names	NUC1
Protein Names	Recommended name: Mitochondrial nuclease EC= 3.1.30.-
Expression Region	1-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.