



# Recombinant Human Rho guanine nucleotide exchange factor 9 (ARHGEF9)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP002057HU   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | O43307   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Homo sapiens (Human)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)  |
| <b>Sequence</b>          | <p>MTLLITGDSI VSAEAVWDHV TMANRELAFK AGDVIKVLDA SNKDWWWGQI<br/> DDEEGWFPAS FVRLWVNQED EVEEGPSDVQ NGHLDPNSDC<br/> LCLGRPLQNR DQMRANVINE IMSTERHYIK HLKDICEGYL KQCRKRRDMF<br/> SDEQLKVIFG NIEDIYRFQM GFVRDLEKQY NNDDPHLSEI GPCFLEHQDG<br/> FWIYSEYCNN HLDACMELSK LMKDSRYQHF FEACRLLQQM IDIAIDGFLL<br/> TPVQKICKYP LQLAELLKYT AQDHSYRYV AAALAVMRNV TQQINERKRR<br/> LENIDKIAQW QASVLDWEGE DILDRSSELI YTGEMAWIYQ PYGRNQQRVF<br/> FLFDHQMVLC KKDILRRDIL YYKGRIDMDK YEVDIEDGR DDDFNVSMKN<br/> AFKLNKETE EIHLFFAKKL EEKIRWLRAF REERKMVQED EKIGFEISEN<br/> QKRQAAMTVR KVPKQKGVNS ARSVPPSYPP PQDPLNHGQY LVPDGLAQSQ<br/> VFEFTEPKRS QSPFWQNFSSR LTPFKK</p> |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | ARHGEF9  |
| <b>Protein Names</b>     | Recommended name: Rho guanine nucleotide exchange factor 9<br>Alternative name(s): Collybistin PEM-2 homolog Rac/Cdc42 guanine nucleotide exchange factor 9  |
| <b>Expression Region</b> | 1-516  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |