



# Recombinant *Saccharomyces cerevisiae* Guanine nucleotide-binding protein alpha-1 subunit (GPA1)

<b>Product Code</b>	CSB-BP362451SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P08539
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GCTVSTQTI GDESDPFLQN KRANDVIEQS LQLEKQRDKN EIKLLLLGAG ESGKSTVLKQ LKLLHQGGFS HQERLQYAQV IWADAIQSMK ILIIQARKLG IQLDCDDPIN NKDLFACKRI LLKAKALDYI NASVAGGSDF LNDYVLKYSE RYETRRRVQS TGRAKAAFDE DGNISNVKSD TDRDAETVTQ NEDADRNNSS RINLQDICKD LNQEGDDQMF VRKTSREIQG QNRRNLIHED IAKAIKQLWN NDKGIKQCFA RSNEFQLEGS AAYYFDNIEK FASPNYVCTD EDILKGRIKT TGITETEFNI GSSKFKVLDA GGQRSERKKW IHCFEGITAV LFVLAMSEYD QMLFEDERVN RMHESIMLFD TLLNSKWFKD TPFILFLNKI DLFEKVKSM PIRKYFPDYQ GRVGDAAEAGL KYFEKIFLSL NKTNKPIYVK RTCATDTQTM KFVLSAVTDL IIQQNLKKIG II
<b>Source</b>	Baculovirus
<b>Target Names</b>	GPA1
<b>Protein Names</b>	Recommended name: Guanine nucleotide-binding protein alpha-1 subunit Alternative name(s): GP1-alpha
<b>Expression Region</b>	2-472
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.