



Recombinant Human Complement C4-A (C4A), partial

Product Code	CSB-MP003949HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C0L4
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	K PRLLLFSPSV VHLGVPLSVG VQLQDVPRGQ VVKGSVFLRN PSRNNVPCSP KVDFTLSSER DFALLSLQVP LKDAKSCGLH QLLRGPEVQL VAHSPWLKDS LSRTTNIQGI NLLFSSRRGH LFLQTDQPIY NPGQRVRYRV FALDQKMRPS TDTITVMVEN SHGLRVRKKE VYMPSSIFQD DFVIPDISEP GTWKISARFS DGLESNSSTQ FEVKKYVLPN FEVKITPGKP YILTVPGHLD EMQLDIQARY IYGKPVQGVA YVRFGLLDED GKKTFFRGLE SQTCLVNGQS HISLSKAEFQ DALEKLNMG I TDLQGLRLYV AAAIIESPGG EMEEAELTSW YFVSSPFLSLD LSKTKRHLVP GAPFLLQALV REMSGSPASG IPVKVSATVS SPGSVPEVQD IQQNTDGSGQ VSIPIIPQT ISELQLSVSA GSPHPAIARL TVAAPPSGGP GFLSIERPDS RPPRVGDTLN LNLRAVGSGA TFSHYYYMIL SRGQIVFMNR EPKRTLTSVS VFVDHHLAPS FYFVAFYYHG DHPVANSLRV DVQAGACEGK LELSVDGAKQ YRNGESVKLH LETDSLALVA LGALDTALYA AGSKSHKPLN MGKVFEAMNS YDLGCGPGGG DSALQVFQAA GLAFSDGDQW TLSRKRLSCP KEKTT
Source	Mammalian cell
Target Names	C4A
Protein Names	Recommended name: Complement C4-AA Alternative name(s): Acidic complement C4 C3 and PZP-like alpha-2-macroglobulin domain-containing protein 2 Cleaved into the following 6 chains: 1. Complement C4 beta chain 2. Complement C4-A alpha cha
Expression Region	20-675
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.