



# Recombinant *Drosophila melanogaster* Homeobox protein abdominal-B (Abd-B)

<b>Product Code</b>	CSB-YP362589DLU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P09087
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila melanogaster</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           MQQHHLQQQQ QQQQQQEQQH LQEQQQHLQQ LHHHAHHHLP            QPLHTTSHHH SAHPHLQQQQ QQQQHAVVAS SPSSVLQQQQ            QQSTPTTHST PTHAVMYEDP PPVPLVAVQQ QHLPAPQQQQ            QLQQQQQQQQ QQLATTPVAG ALSPAQTPTG PSAQQQQHLT            SPHHQQLPQQ QTPNSVASGA SSNLQQQQQQ QNAAVAPGQT            QIVAPTTASV SPSSVSSQKE DINMSIQLAP LHIPAIRAGP GFETDTSAAV            KRHTAHWAYN DEGFNQHYGS GYYDRKHMFA YYPETQFPV            GQYWGPNYRP DQTTSAAAAA AYMNEAERHV SAAARQSVEG            TSTSSYEPPT YSSPGGLRGY PSENYSSSGA SGGLSVGAVG PCTPNPGLHE            WTGQVSVRKK RKPYSKFQTL ELEKEFLFNA YVSKQKRWEL ARNLQLTERQ            VKIWFQNRMM KNKKNSQRQA NQQNNNNNSS SNHNHAQATQ            QHHSQHHLNL SLNMGHHAAK MHQ         </p>
<b>Source</b>	Yeast
<b>Target Names</b>	Abd-B
<b>Protein Names</b>	Recommended name: Homeobox protein abdominal-B Alternative name(s): Infraabdominal 7 Short name= IAB-7 P3 PH189
<b>Expression Region</b>	1-493
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.