



Recombinant Enterobacteria phage 933W Shiga-like toxin 2 subunit A (stxA2)

Product Code	CSB-EP362675ECC-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09385
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Enterobacteria phage 933W (Bacteriophage 933W)
Purity	>85% (SDS-PAGE)
Sequence	REFTIDFS TQQSYVSSLN SIRTEISTPL EHSQGTTSV SVINHTPPGS YFAVDIRGLD VYQARFDHLR LIIEQNNLYV AGFVNTATNT FYRFSDFTHI SVPGVTTVSM TTDSSYTTLQ RVAALERSGM QISRHSLVSS YLALMEFSGN TMTRDASRAV LRFVTVTAEA LRFRQIQREF RQALSETAPV YTMTPGDVDL TLNWGRISNV LPEYRGEDGV RVGRISFNNI SAILGTVAVI LNCHHQGARS VRAVNEESQP ECQITGDRPV IKINNTLWES NTAAAF LN RK SQFLYTTGK
Source	E.coli
Target Names	stxA2
Protein Names	Recommended name: Shiga-like toxin 2 subunit A Short name= SLT-2 A subunit Short name= SLT-2a Short name= SLT-IIa EC= 3.2.2.22Alternative name(s): Verocytotoxin 2 subunit A Verotoxin 2 subunit A rRNA N-glycosid
Expression Region	23-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.