



# Recombinant Human Ribosomal RNA-processing protein 8 (RRP8)

<b>Product Code</b>	CSB-EP020531HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O43159
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFEEPEWAEA APVAAGLGPV ISRPPPAASS QNKGSKRRQL LATLRALEAA SLSQHPPSLC ISDSEEEEE RKKKCPKKAS FASASAEVGK KGKKKCQKQG PPCSDSEEEV ERKKKCHKQA LVGSDSAEDE KRKRKCQKHA PINSAQHLDN VDQTGPKAWK GSTTNDPPKQ SPGSTSPKPP HTLSRKQWRN RQKNKRCKN KFQPPQVPDQ APAEAPTEKT EVSPVPRDTS HEARAGALRA RMAQRDLGAR FRYLNEQLYS GPSSAAQRLF QEDPEAFLLY HRGFQSQVKK WPLQPVDRIA RDLRQRPASL VVADFGCGDC RLASSIRNPV HCFDLASLDP RVTVCDMAQV PLEDESVDVA VFCLSLMGTN IRDFLEEANR VLKPGGLLKV AEVSSRFEDV RTFLRAVTKL GFKIVSKDLT NSHFFLDFDQ KTGPPPLVGPK AQLSGLQLQP CLYKRR
<b>Source</b>	E.coli
<b>Target Names</b>	RRP8
<b>Protein Names</b>	Recommended name: Ribosomal RNA-processing protein 8 EC= 2.1.1.- Alternative name(s): Cerebral protein 1 Nucleomethylin
<b>Expression Region</b>	1-456
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.