



Recombinant Zea mays DNA-directed RNA polymerase subunit alpha (rpoA)

Product Code	CSB-MP362730ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09562
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	MVREEITGST QTLEWKCVES RVDSKRLYYG RFILSPLRKG QADTVGIALR RALLGEIEGT CITRAKFGNV PHEYSTIVGI EESIQEILLN LKEIVLRSNL YGVRDASICV KGPRYITAQD IILPPSVEIV DTTQPIANLR EPVDFCIELQ IKRDRAYHTE LRKNSQDGSY PIDAVFMPVR NVNYSIFSCG NGNEKHEILF LEIWTNGSLT PKEALYEASR NLIDLFLPFL HTEEEGTSFE ENKNRLTPPL LTFQKRFTNL KKNKKGIPLN CIFIQLELP SRTYNCLKRA NIHTLLDLLS KTEEDLMQIN SFRMEDGKLI WDTLEKHLPI DLPKNKFSL
Source	Mammalian cell
Target Names	rpoA
Protein Names	Recommended name: DNA-directed RNA polymerase subunit alpha Short name= PEP EC= 2.7.7.6 Alternative name(s): Plastid-encoded RNA polymerase subunit alpha Short name= RNA polymerase subunit alpha
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.