



Recombinant Sphingomyelinase C (sph)

Product Code	CSB-MP362741BQJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09599
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus
Purity	≥85% (SDS-PAGE)
Sequence	DTS TDQNNTLKVM THNVYMLSTN LYPNWGQSQR ADLIGAADYI KNQDVVILNE VFDNSASDRL LGNLKKEYPN QTAVLGRSNG NEWDKTLGSY SSSTPEDGGV AIVSKWPIVE KIQYVFAKGC GPDNLSNKGK VYTKIKKNDR FVHVIGTHLQ AEDSMCGKTS PASVRTNQLK EIQDFIKNKN IPNDEYVLFQ GDMNVNKINA ENNSDSEYAS MFKTLHASIP SYTGHTATWD ATTNIAKYN FPDSPAAYLD YIIASKDHAN PSFIENKVLQ PKSPQWTVTS WLKKYTYDDY SDHYPVAATI SMK
Source	Mammalian cell
Target Names	sph
Protein Names	Recommended name: Sphingomyelinase C Short name= SMase EC= 3.1.4.12 Alternative name(s): Cereolysin B SMPLC Sphingomyelin phosphodiesterase
Expression Region	28-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.