



Recombinant Enterobacteria phage SP Readthrough protein A1

Product Code	CSB-MP362757EDU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09677
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage SP (Bacteriophage SP)
Purity	≥85% (SDS-PAGE)
Sequence	MAKLNQVTLK KIGKNGDQTL TLTPRGVNPT NGVASLSEAG AVPALEKRVT VSVAQPSRNR KNFKVQIKLQ NPTACTRDAC DPSVTRSAFA DVTLSTSYS TDEERALIRT ELAALLADPL IVDAIDNLNP AYWAALLVAS SGGGDNPSPD DVPVVPDVKP PDGTGRYKCP FACYRLGSIY EVGKEGSPDI YERGDEVSVT FDYALEDFLG NTNWRNWDQR LSDYDIANRR RCRGNGYIDL DATAMQSDDF VLSGRYGVRK VKFPGAFGSI KYLLNIQGDA WLDLSEVTAY RSYGMVIGFW TDSKSPQLPT DFTQFNSANC PVQTVIIIPS L
Source	Mammalian cell
Protein Names	Recommended name: Readthrough protein A1 Cleaved into the following chain: 1. Coat protein
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.