



Recombinant 50S ribosomal protein L1 (rpIA)

Product Code	CSB-EP362786SYN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09764
Product Type	Recombinant Protein
Immunogen Species	Serratia marcescens
Purity	>85% (SDS-PAGE)
Sequence	AKLTKRMRV IRDKVDATKQ YDITEAIALL KELATAKFVE SVDVAVNLGI DARKSDQNVR GATVLPHTG RSVRVAVFTQ GANAEAAKAA GAELVGMEDL AEQIKKGEMN FDVVIASPDA MRVVGQLGQI SGPRGLMPNP KVGTVTPNVA EAVKNAKAGQ VRYRNDKNGI IHTTIGKVDF DADKLKENLE ALLVALKKAK PSQAKGMYIK KVSLSTTMGA GVAIDQSGLS AAAN
Source	E.coli
Target Names	rpIA
Protein Names	Recommended name: 50S ribosomal protein L1
Expression Region	2-234
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.