



Recombinant Neisseria meningitidis serogroup B Outer membrane protein class 4 (rmpM)

Product Code	CSB-BP362992NKG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A0V3
Product Type	Recombinant Protein
Immunogen Species	Neisseria meningitidis serogroup B (strain MC58)
Purity	>85% (SDS-PAGE)
Sequence	GEASVQGY TVSGQSNEIV RNNYGECWKN AYFDKASQGR VECGDAVAAP EPEPEPEPAP APVVVVEQAP QYVDETISLS AKTLFGFDKD SLRAEAQDNL KVLAQRLSRT NVQSVRVEGH TDFMGSDKYN QALSERRAYV VANNLVSNGV PVSRI SAVGL GESQAQMTQV CEAEVAKLGA KVS KAKKREA LIACIEPDRR VDVKIRSIVT RQV VPAHNHH QH
Source	Baculovirus
Target Names	rmpM
Protein Names	Recommended name: Outer membrane protein class 4
Expression Region	23-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.