



Recombinant Surface presentation of antigens protein spaN (spaN)

Product Code	CSB-EP363108SZB-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A1K5
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	ALDNINLNF SSDKQIEKCE KLSSIDNIDS LVLKKKRKVE IPEYSLIASN YFTIDKHFEH KHDKGEIYSG IKNAFELRNE RATYSDIPES MAIKENILIP DQDIKAREKI NIGDMRGIFS YNKSGNADKN FERSHTSSVN PDNLLESDNR NGQIGLKNHS LSIDKNIADI ISLLNGSVAK SFELPVMNKN TADITPSMSL QEKSIENDK NVFQKNSEMT YHFKQWGAGH SVSISVESGS FVLKPSDQFV GNKLDLILKQ DAEGNYRFDS SQHNKGNKNN STGYNEQSEE EC
Source	E.coli
Target Names	spaN
Protein Names	Recommended name: Surface presentation of antigens protein spaN Alternative name(s): Spa32 protein
Expression Region	2-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.