



Recombinant Mouse Guanylate cyclase soluble subunit beta-1 (Gucy1b3)

Product Code	CSB-EP010052MO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O54865
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MYGFVNHLE LLVIRNYGPE VWEDIKKEAQ LDEEGQFLVR IYDDSKTYD LVAASKVLN LNAGEILQMF GKMFFVFCQE SGYDTILRVL GSNVREFLQN LDALHDHLAT IYPGMRAPSF RCTDAEKGGK LILHYYSERE GLQDIVIGII KTVAQQIHGT EIDMKVIQQR NEECDHTQFL IEEKESKEED FYEDLDRFEE NGTQESRISP YTFCKAFPFFH IIFDRNLVVT QCGNAIYRVL PQLQPGNCSL LSVFSLVRPH IDISFHGILS HINTVFLRS KEGLLDVEKL ECEDELTAE ISCLRLKGQM IYLPEADSIL FLCSPSVMNL DDLTRRGLYL SDIPLHDATR DLVLLGEQFR EYKLTQELE ILTDRQLTL RALEDEKKKT DTLLYSVLPP SVANELRHKR PVPARYDNV TILFSGIVGF NAFCSKHASG EGAMKIVNLL NDLYTRFDL TDSRKNPFVY KVETVGDYKYM TVSGLPEPCI HHARSICHLA LDMMEIAGQV QVDGESVQIT IGIHTGEVVT GVIGQRMPRY CLFGNTVNLT SRTETTGEKG KINVSEYTYR CLMSPENSDP LFHLEHRGPV SMKGKKEPMQ VWFLSRKNTG TEETNEEDEN</p>
Source	E.coli
Target Names	Gucy1b1
Protein Names	Recommended name: Guanylate cyclase soluble subunit beta-1 Short name= GCS-beta-1 EC= 4.6.1.2 Alternative name(s): Guanylate cyclase soluble subunit beta-3 Short name= GCS-beta-3 Soluble guanylate cyclase small subunit
Expression Region	1-620
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.