



Recombinant Mouse Nuclear receptor subfamily 1 group I member 2 (Nr1i2)

Product Code	CSB-MP016048MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O54915
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MRPEESWSRV GLVQCEEADS ALEEPINVEE EDGGLQICRV CGDKANGYHF NVMTCEGCKG FFRRAMKRNV RLRCPFRKGT CEITRKTRRQ CQACRLRKCL ESGMKKEMIM SDAAVEQRRRA LIKRKKREKI EAPPPGGQGL TEEQQALIQE LMDAQMOTFD TTFSHFKDFR LPAVFHSGCE LPEFLQASLL EDPATWSQIM KDRVPMKISL QLRGEDGSIW NYQPPSKSDG KEIIPLLPHL ADVSTYMFKG VINFAKVISY FRDLPIEDQI SLLKGATFEM CILRFNTMFD TETGTWECGR LAYCFEDPNG GFQKLLLDPL MKFHCMLKKL QLHKEEYVLM QAISLFSPDR PGVVQRSVVD QLQERFALTL KAYIECSRPY PAHRFLFLKI MAVLTELRSI NAQQTQQLLR IQDSHPFATP LMQELFSSTD G
Source	Mammalian cell
Target Names	Nr1i2
Protein Names	Recommended name: Nuclear receptor subfamily 1 group I member 2 Alternative name(s): Orphan nuclear receptor PXR Pregnane X receptor
Expression Region	1-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene product belongs to the nuclear receptor superfamily, members of which are transcription factors characterized by a ligand-binding domain and a DNA-binding domain. The encoded protein is a transcriptional regulator of the cytochrome P450 gene CYP3A4, binding to the response element of the CYP3A4 promoter as a heterodimer with the 9-cis retinoic acid receptor RXR. It is activated by a range of compounds that induce CYP3A4, including dexamethasone and rifampicin. Several alternatively spliced transcripts encoding different isoforms, some of which use non-AUG (CUG) translation initiation codon, have been described for this gene. Additional transcript variants exist, however, they have not been fully characterized.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.