



Recombinant Salmonella typhimurium Phosphate acyltransferase (plsX)

Product Code	CSB-MP363201SXB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A259
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Purity	>85% (SDS-PAGE)
Sequence	MTRLTLALDV MGGDFGPSVT VPAALQALNA NSQLTLLLVG NPDIITPLLA KADFEQRSRL QIIPAQSVIA SDARPSQAIR ASRGTSMRVA LELVKEGRAE ACVSAGNTGA LMGLAKLLLK PLEGIERPAL VTVLPHQQKG KTVVLDLGAN VDCDSTMLVQ FAVMGAVLAE EVVGIKNPRV ALLNIGEEET KGLDSIREAS LMLKTVPTIN YIGYLEANEL LTGKTDVLVC DGFTGNVTLK TMEGVVVMFL SLLKSQGEGK KRSWWLLLLK RWLQKSLTRR FSHLNPQYN GACLLGLRGT VIKSHGAANQ RAFAVAIEQA VQAVQRQVPQ RIAARLESVY PAGFEPLDDG KGVNLRHR
Source	Mammalian cell
Target Names	plsX
Protein Names	Recommended name: Phosphate acyltransferase EC= 2.3.1.n2 Alternative name(s): Acyl-ACP phosphotransacylase Acyl-[acyl-carrier-protein]--phosphate acyltransferase Phosphate-acyl-ACP acyltransferase
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.