



Recombinant Human Neuropilin-2 (NRP2)

Product Code	CSB-BP016092HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O60462
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	RGQPDPPCGGRLNSKDAGYITSPGYDYPHQNCEWIVYAPEPNQKIVLNF NPHFEIEKHDKYDFIEIRDGDSEADLLGKHCGNIAPPTIISGSMYIKFTSDY ARQGAGFSLRYEIFKTGSEDCSKNFTSPNGTIESPGFPEKYPHNLDCTFTILAK PKMEIILQFLIFDLEHDPLQVGECDKYDWLDIWDGIPHVGPLIGKYCGTKTPS ELRSSTGILSLTFHTDMAVAKDGFSAARYLVHQEPLNFQCNVPLGMESGRIA NEQISASSTYS DGRWTPQQSRLHGDDNGWTPNLDSNKEYLQVDLRLTMLTA IATQGAISRETQNGYYVKSYLEVSTNGEDWMVYRHGKNHKVVFQANNDATEV VLNKLHAPLLTRFVRIRPQTWHSIALRLELFGCRVTDAPCSNMLGMLSLGIAD SQISASSTQEYLWSPSAARLVSSRSGWFPRIQAQPGEEWLQVDLGTGPKTVK GVIIQGARGGDSITAVEARAFVRKFKVSYSLNGKDWEYIQDPRTQQPKVGC WRPL
Source	Baculovirus
Target Names	NRP2
Protein Names	Recommended name: Neuropilin-2 Alternative name(s): Vascular endothelial cell growth factor 165 receptor 2
Expression Region	21-555aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein of Isoform s9
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.