



Recombinant Mycobacterium tuberculosis MPT51/MPB51 antigen (mpt51)

Product Code	CSB-MP363615MVZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P9WQN6
Product Type	Recombinant Protein
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Purity	>85% (SDS-PAGE)
Sequence	AEPT AKAAPYENLM VPSPSMGRDI PVAFLAGGPH AVYLLDAFNA GPDVSNWVTA GNAMNTLAGK GISVVAPAGG AYSMYTNWEQ DGSKQWDTFL SAELPDWLAA NRGLAPGGHA AVGAAQGGYG AMALAAFHPD RFGFAGSMSG FLYPSNTTTN GAIAAGMQQF GGVDTNGMWG APQLGRWKWH DPWVHASLLA QNNTRVWVWS PTNPGASDPA AMIGQAAEAM GNSRMFYNQY RSVGGHNGHF DFPASGDNGW GSWAPQLGAM SGDIVGAIR
Source	Mammalian cell
Target Names	mpt51
Protein Names	Recommended name: MPT51/MPB51 antigen
Expression Region	27-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.