




# Recombinant Mycobacterium tuberculosis ESAT-6-like protein esxB (esxB)

<b>Product Code</b>	CSB-YP363658MVZ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P9WNK4
<b>Alias</b>	10 kDa culture filtrate antigen CFP-10 Secreted antigenic protein MTSA-10
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	AEMKTDAAATLAQEAGNFERISGDLKTQIDQVESTAGSLQGQWRGAAGTAAQA AVVRFQEAANKQKQELDEISTNIRQAGVQYSRADEEQQQALSSQMGE
<b>Lead Time</b>	3-7 business days
<b>Research Area</b>	Others
<b>Source</b>	Yeast
<b>Gene Names</b>	esxB
<b>Protein Names</b>	Recommended name: ESAT-6-like protein esxB Alternative name(s): 10 kDa culture filtrate antigen CFP-10 Secreted antigenic protein MTSA-10
<b>Expression Region</b>	2-100aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	12.7kDa
<b>Protein Description</b>	Full Length of Mature Protein

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.