




Recombinant Mycobacterium tuberculosis ESAT-6-like protein esxB (esxB)

Product Code	CSB-YP363658MVZ
Abbreviation	Recombinant Mycobacterium tuberculosis esxB protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P9WNK4
Alias	10 kDa culture filtrate antigen CFP-10 Secreted antigenic protein MTSA-10
Product Type	Recombinant Protein
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	AEMKTDAAATLAQEAGNFERISGDLKTQIDQVESTAGSLQGQWRGAAGTAAQA AVVRFQEAANKKQKQELDEISTNIRQAGVQYSRADEEQQQALSSQMGE
Research Area	Others
Source	Yeast
Target Names	esxB
Protein Names	Recommended name: ESAT-6-like protein esxB Alternative name(s): 10 kDa culture filtrate antigen CFP-10 Secreted antigenic protein MTSA-10
Expression Region	2-100aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	12.7kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.