



Recombinant Acetate kinase (ackA)

Product Code	CSB-BP363835EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A6A4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MSSKLVLVLN CGSSSLKFAI IDAVNGEEYL SGLAECFHLP EARIKWKMDG NKQEAALGAG AAHSEALNFI VNTILAQKPE LSAQLTAIGH RIVHGGEKYT SSVIDESVI QGIKDAASFA PLHNPALHIG IEEALKSFPQ LKDKNVAVFD TAFHQTMPEE SYLYALPYNL YKEHGIRRYG AHGTSHFYVT QEAAKMLNKP VEELNIITCH LGNGGSVSAI RNGKCVDTSM GLTPLEGLVM GTRSGDIDPA IIFHLHDTLG MSVDAINKLL TKESGLLGLT EVTSDCRYVE DNYATKEDAK RAMDVYCHRL AKYIGAYTAL MDGRLDAVVF TGGIGENAAM VRELSLGKLG VLGFEVDHER NLAARFGKSG FINKEGTRPA VVIPTNEELV IAQDASRLTA
Source	Baculovirus
Target Names	ackA
Protein Names	Recommended name: Acetate kinase EC= 2.7.2.1 Alternative name(s): Acetokinase
Expression Region	1-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.