



# Recombinant Peptide chain release factor 3 (prfC)

<b>Product Code</b>	CSB-MP364010EGX
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0A7I5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TLSPYLQEV AKRRTFAIIS HPDAGKTTIT EKVLLFGQAI QTAGTVKGRG SNQHAKSDWM EMEKQRGISI TTSVMQFPYH DCLVNLLDTP GHEDFSEDY RTLTAVDCCL MVIDAAKGVE DRTRKLMEVT RLRDTPILTF MNKLDRDIRD PMELLDEVEN ELKIGCAPIT WPIGCGKLFK GVYHLYKDET YLYQSGKGHT IQEVRIKGL NNPDLDAAVG EDLAQQLRDE LELVKGASNE FDKELFLAGE ITPVFFGTAL GNFGVDHMLD GLVEWAPAPM PRQTDTRTVE ASEDKFTGFV FKIQANMDPK HRDRVAFMRV VSGKYEKGMK LRQVRTAKDV VISDALTFMA GDRSHVEEAY PGDILGLHNH GTIQIGDTFT QGEMMKFTGI PNFAPELFRR IRLKDPLKQK QLLKGLVQLS EEGAVQVFRP ISNNDLIVGA VGVLLQFDVVV ARLKSEYNVE AVYESVNVAT ARWVECADAK KFEFEKCRKNE SQLALDGGDN LAYIATSMVN LRLAQERYPD VQFHQTREH
<b>Source</b>	Mammalian cell
<b>Target Names</b>	prfC
<b>Protein Names</b>	Recommended name: Peptide chain release factor 3 Short name= RF-3
<b>Expression Region</b>	2-529
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.