



Recombinant S-adenosylmethionine synthase (metK)

Product Code	CSB-BP364097EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A818
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	AKHLFTSES VSEGHPDKIA DQISDAVLDA ILEQDPKARV ACETYVKTGM VLVGGEITTS AWVDIEEITR NTVREIGYVH SDMGFDANSC AVLSAIGKQS PDINQGV DRA DPLEQGAGDQ GLMFGYATNE TDVLMPAPIT YAHRLVQRQA EVRKNGTLPW LRPDAKSQVT FQYDDGKIVG IDAVVLSTQH SEEIDQKSLQ EAVMEEIIPK ILPAEWL TSA TKFFINPTGR FVIGGPMGDC GLTGRKIIVD TYGGMARHGG GAFSGKDPSK VDRSAAYAAR YVAKNIVAAG LADRCEIQVS YAIGVAEPTS IMVETFGTEK VPSEQLTLLV REFFDLRPYG LIQMLDLLHP IYKETAAYGH FGREHFPWEK TDKAQLLRDA AGLK
Source	Baculovirus
Target Names	metK
Protein Names	Recommended name: S-adenosylmethionine synthase Short name= AdoMet synthase EC= 2.5.1.6 Alternative name(s): MAT Methionine adenosyltransferase
Expression Region	2-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.