



# Recombinant Human Suppressor of cytokine signaling 5 (SOCS5)

<b>Product Code</b>	CSB-MP022394HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O75159
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDKVGKMWNN FKYRCQNLFG HEGGSRSENV DMNSNRCLSV KEKNISIGDS TPQQQSSPLR ENIALQLGLS PSKNSSRRNQ NCATEIPQIV EISIEKDNDSD CVTPGTRLAR RDSYSRHAPW GGKKKHSCST KTQSSLDADK KFGRTSGLQ RRERRYGVSS VHDMDSVSSR TVGSRSLRQR LQDTVGLCFP MRTYSKQSKP LFSNKRKIHL SELMLEKCPF PAGSDLAQKW HLIKQHTAPV SPHSTFFDTF DPSLVSTEDE EDRLRERRRL SIEEGVDPPP NAQIHTFEAT AQVNPLYKLG PKLAPGMTEI SGDSSAIPQA NCDSEEDTTT LCLQSRRQKQ RQISGDSHTH VSRQGAWKVH TQIDYIHCLV PDLLQITGNP CYWGVMDRYE AEALLEGKPE GTFLLRDSAQ EDYLFSVSFR RYNRSLHARI EQWNHNFSFD AHDPCVFHSS TVTGLLEHYK DPSSCMFFEP LLTISLNRTF PFSLQYICRA VICRCTTYDG IDGLPLPSML QDFLKEYHYK QKVRVRWLER EPVKAK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SOCS5
<b>Protein Names</b>	Recommended name: Suppressor of cytokine signaling 5 Short name= SOCS-5 Alternative name(s): Cytokine-inducible SH2 protein 6 Short name= CIS-6 Cytokine-inducible SH2-containing protein 5
<b>Expression Region</b>	1-536
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This protein contains a SH2 domain and a SOCS BOX domain. The protein thus belongs to the suppressor of cytokine signaling (SOCS) family, also known as STAT-induced STAT inhibitor (SSI) protein family. SOCS family members are known to be cytokine-inducible negative regulators of cytokine signaling. The specific function of this protein has not yet been determined. Two alternatively spliced transcript variants encoding an identical protein have been reported.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.