



# Recombinant Serine--tRNA ligase (serS)

<b>Product Code</b>	CSB-EP364184EOD
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0A8L3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLDPNLLRNE PDAVAEKLAR RGFKLDVDKL GALEERRKVL QVKTENLQAE RNSRSKSIQ AKARGEDIEP LRLEVNKLGE ELDAAKAELD ALQAEIRDIA LTIPNLPAD VPVGKDENDN VEVSRWGTTPR EFD FEVRDHV TLGEMHSGLD FAAAVKLTGS RFVVMKGQIA RMHRALSQFM LDLHTEQHG Y SENYVPYLVN QDTLYGTGQL PKFAGDLFHT RPLEEEADTS NYALIPTAEV PLTNLVRGEI IDEDDLPIKM TAHTPCFRSE AGSYGRDTRG LIRMHQFDKV EMVQIVRPED SMAALEEMTG HAEKVLQLLG LPYRKIILCT GDMGFGACKT YDLEVVWIPAQ NTYREISSCS NVWDFQARRM QARCRSKSDK KTRLVHTLNG SGLAVGRTL V AVMENYQQAD GRIEVPEVLR PYMNGLEYIG
<b>Source</b>	E.coli
<b>Target Names</b>	serS
<b>Protein Names</b>	Recommended name: Serine--tRNA ligase EC= 6.1.1.11 Alternative name(s): Seryl-tRNA synthetase Short name= SerRS Seryl-tRNA(Ser/Sec) synthetase
<b>Expression Region</b>	1-430
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.