



Recombinant KHG/KDPG aldolase (eda)

Product Code	CSB-MP364279EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A956
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MKNWK TSAES ILTTGPVVPV IVVKLEHAV PMAKALVAGG VRVLEVT LRT ECAVD AIRAI AKEVPEAIVG AGTVLNPQQL AEVTEAGA QF AISPLGTEPL LKAATEGTIP LIPGISTVSE LMLGMDYGLK EFKFFPAEAN GGVKALQAIA GPF SQVRFCP TGGISPANYR DY LALKSVLC IGGSWLV PAD ALEAGDYDRI TKLAREAVEG AKL
Source	Mammalian cell
Target Names	eda
Protein Names	Recommended name: KHG/KDPG aldolase Including the following 2 domains: 4-hydroxy-2-oxoglutarate aldolase EC= 4.1.3.16 Alternative name(s): 2-keto-4-hydroxyglutarate aldolase Short name= KHG-aldolase 2-dehydro-3-deoxy-phos
Expression Region	1-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.