



# Recombinant Escherichia coli Anaerobic glycerol-3-phosphate dehydrogenase subunit A (glpA)

<b>Product Code</b>	CSB-EP364308ENV-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0A9C0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKTRDSQSSD VIIIIGGGATG AGIARDCALR GLRVILVERH DIATGATGRN HGLLHSGARY AVTDAESARE CISENQILKR IARHCVEPTN GLFITLPEDD LSFQATFIRA CEEAGISAEA IDPQQARIIE PAVNPALIGA VKVPDGTVDP FRLTAAANMLD AKEHGAVILT AHEVTGLIRE GATVCGVRVR NHLTGETQAL HAPVVVNAAG IWGQHIAEYA DLRIRMFPK GLLIMDHRI NQHVINRCRK PSDADILVPG DTISLIGTTS LRIDYNEIDD NRVTAAEVDI LLREGEKLAP VMAKTRILRA YSGVRPLVAS DDDPSGRNVS RGIVLLDHAE RDGLDGFITI TGGKLMTYRL MAEWATDAVC RKLGNTRPCT TADLALPGSQ EPAEVTLRKV ISLPAPLRGS AVYRHGDRTP AWLSEGR LHR SLVCECEAVT AGEVQYAVEN LNVNSLLDLR RRTRVGMGTC QGELCACRAA GLLQRFNVTT SAQSIEQLST FLNERWKG VQ PIAWGDALRE SEFTRWVYQG LCGLEKEQKD AL
<b>Source</b>	E.coli
<b>Target Names</b>	glpA
<b>Protein Names</b>	Recommended name: Anaerobic glycerol-3-phosphate dehydrogenase subunit A Short name= G-3-P dehydrogenase EC= 1.1.5.3
<b>Expression Region</b>	1-542
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.